



BIPOLAR  
FORCEPS

**Kali**®

Single Use  
Bipolar Forceps





## Kali®

### Single use bipolar forceps

#### Technical features

- Optimal tissue cauterization for a safer surgery
- Grasping and handling of tissue without changing instrument
- Fine dissection allowed by specific shape of jaws
- Operational performance at all times
- 360° rotation of jaws
- Preconnected to a single use bipolar cable (3-meter length – Coax 4/8 connection)
- Sterile, ready to use

Product reference	Designation	Length	Diameter	Units per box
BP5330	«Kali®» Bipolar forceps	330 mm	5 mm	5

MD Class IIb.  
Sterilised with ETO  
See instructions for use



#### Manufactured by:

Péters Surgical  
42 rue Benoit Frachon  
93013 Bobigny Cedex - FRANCE

Tel: +33 (0)1 48 10 62 62  
Fax: +33 (0)1 48 91 22 99  
commercial@peters-surgical.com  
www.peters-surgical.com